

Direct-add bacteria can be added directly pouring the acclimated bag into the wine. Alternatively, they can be propagated with Reactivateur to increase population and efficiency:

1. Draw a small portion of the wine to be inoculated.
2. Use 100 liters (26 Gal) for 250 hL (6600 Gal).
3. Add 60 ppm of Fermoplus Malolactique.
4. Adjust pH to 3.5 – 4 and inoculate with the malolactic bacteria of your choice.
5. Maintain a constant temperature of 24°C/75°F for 24 hrs.
6. The next day that portion of wine will have a much more aggressive population to quickly start the ML process in the rest of the tank.